



Compatibility Chart of INNO Implant System ®

| Cowellmedi | | | Astra | | | | Straumann | | | | Nobel Biocare | | | | 3i | | | | | |
|-----------------------------|----------|---------------------|---------|----------|---------------------|---------------|------------------|-----------------|----------|---------------|---------------|---------------|---------------------|---------------|---------|-------------------------|---------------------|---------------|---------|-----|
| Product | Diameter | Internal Connection | Product | Diameter | Internal Connection | Compatibility | Product | Diameter | Platform | Compatibility | Product | Diameter | Internal Connection | compatibility | Product | Diameter | Internal Connection | Compatibility | | |
| INNO Sub. Narrow | 3.3 | 2.1 Hex | | | | | | | | | | | | | | | | | | |
| INNO Sub. / INNO Sub. Short | 3.5 | 2.5 Hex | Yellow | 3.0 | | N | St. Taper (Int.) | 3.3 | 3.5 | N | Nobel Active | | | N | Prevail | | | N | | |
| | 4.0 | 2.5 Hex | Aqua | 3.5 | | N | | 3.3 | 4.8 | Y | | | | N | | | N | | | |
| | 4.5 | 2.5 Hex | | 4.0 | | N | | 4.1 | 4.8 | Y | | | | N | | | N | | | |
| | 5.0 | 2.5 Hex | Lilac | 4.5 | | Y | | 4.8 | 4.8 | Y | | | | N | | | N | | | |
| | 6.0 | 2.5 Hex | | 5.0 | | Y | | 4.8 | 6.5 | Y | | | | N | | | N | | | |
| INNO Int. | 3.5 | 3.1 Octa | | | | | | Only Solid Type | | | | Nobal Replace | | | N | Trans Gingival | | | N | |
| | 4.0 | 3.1 Octa | | | | | | | | N | | | | N | | | N | | | |
| | 4.5 | 3.1 Octa | | | | | | | | N | | | | N | | | N | | | |
| | 5.0 | 3.1 Octa | | | | | | | | N | | | | N | | | N | | | |
| | 6.0 | 3.1 Octa | | | | | | | | N | | | | N | | | N | | | |
| INNO Ext. | 3.5 | 2.7 Hex | | | | | | | | | | Mk III | 3.3 | 2.4 Hex | N | Nanotite (Ext. Implant) | 3.3 | 2.4 Hex | N | |
| | 4.0 | 2.7 Hex | | | | | | | | | | | 3.8 | 2.7 Hex | Y | | | 3.8 | 2.7 Hex | 0.0 |
| | 4.5 | 2.7 Hex | | | | | | | | | | | 4.0 | 2.7 Hex | Y | | | 4.0 | 2.7 Hex | 0.0 |
| | 5.0 | 3.4 Hex | | | | | | | | | | | 5.0 | 3.4 Hex | Y | | | 5.0 | 3.4 Hex | 0.0 |
| | 6.0 | 3.4 Hex | | | | | | | | | | | 5.0 | 3.4 Hex | Y | | | 6.0 | 2.7 Hex | x |

| Cowellmedi | | | Camlog | | | | Hiossen | | | | Dio Implant | | | | MIS | | | | |
|-----------------------------|----------|---------------------|--------------|----------|---------------------|---------------|----------------------------|----------|---------------------|-----------------------|------------------------------------|----------|---------------------|---------------|-----------|----------|---------------------|---------------|---|
| Product | Diameter | Internal Connection | Product | Diameter | Internal Connection | Compatibility | Product | Diameter | Internal Connection | Compatibility | Product | Diameter | Internal Connection | Compatibility | Product | Diameter | Internal Connection | Compatibility | |
| INNO Sub. Narrow | 3.3 | 2.1 Hex | | | | | GSII, GSIII TSII, TSIII | 3.5 | 2.1 Hex | x | | | | | | | | | |
| INNO Sub. / INNO Sub. Short | 3.5 | 2.5 Hex | Tube-In-Tube | | | N | | 4 | 2.5 Hex | Y | UF Impalnt UF II Implant | 3.8 | 2.5 Hex | Y | Seven, M4 | | | | N |
| | 4.0 | 2.5 Hex | | | | N | | 4.5 | 2.5 Hex | Y | | 4.0 | 2.5 Hex | Y | | | | | N |
| | 4.5 | 2.5 Hex | | | | N | | 5 | 2.5 Hex | Y | | 4.5 | 2.5 Hex | Y | | | | | N |
| | 5.0 | 2.5 Hex | | | | N | | 6 | 2.5 Hex | Y | | 5.0 | 2.5 Hex | Y | | | | | N |
| | 6.0 | 2.5 Hex | | | | N | | 7 | 2.5 Hex | Y | | 5.5 | 2.5 Hex | Y | | | | | N |
| INNO Int. | 3.5 | 3.1 Octa | | | | | SS II | 4 | 2.9 Octa | x | Steady Implant system (1 stage) | 3.75 | 2.9 Octa | N | | | | | |
| | 4.0 | 3.1 Octa | | | | 4.5 | | 2.9 Octa | x | 4.0 | | 2.9 Octa | N | | | | | | |
| | 4.5 | 3.1 Octa | | | | 5 | | 2.9 Octa | x | 5.0 | | 2.9 Octa | N | | | | | | |
| | 5.0 | 3.1 Octa | | | | | | | | 5.0 (platform 6.0) | | 2.9 Octa | N | | | | | | |
| | 6.0 | 3.1 Octa | | | | | | | | | | | | | | | | | |
| INNO Ext. | 3.5 | 2.7 Hex | | | | | US | 3.5 | 2.4 Hex | x | FTN /FSN Implant | 3.3 | 2.4 Hex | N | | | | | |
| | 4.0 | 2.7 Hex | | | | 4 | | 2.7 Hex | Y | 3.8 | | 2.7 Hex | Y | | | | | | |
| | 4.5 | 2.7 Hex | | | | 4.5 | | 2.7 Hex | Y | 4.0 | | 2.7 Hex | Y | | | | | | |
| | 5.0 | 3.4 Hex | | | | 5 | | 3.4 Hex | Y | 5.0 | | 3.4 Hex | Y | | | | | | |
| | 6.0 | 3.4 Hex | | | | 5(3i) | | 2.7 Hex | x | 5.5 / 6.0 | | 3.4 Hex | Y | | | | | | |
| | | | | | | | | | | | 5.0(3i) | 2.7 Hex | N | | | | | | |

| Cowellmedi | | | Megagen | | | | Dentis | | | |
|-----------------------------|----------|---------------------|------------------------------|----------|---------------------|---------------|---------|----------|---------------------|---------------|
| Product | Diameter | Internal Connection | Compatibility | Diameter | Internal Connection | Compatibility | Product | Diameter | Internal Connection | Compatibility |
| INNO Sub. Narrow | 3.3 | 2.1 Hex | | | | | | | | |
| INNO Sub. / INNO Sub. Short | 3.7 | 2.5 Hex | Any one | 3.0/3.25 | 2.5 Hex | Y | s-clean | 3.7 | 2.5 Hex | Y |
| | 4.0 | 2.5 Hex | | 4.0 | 2.5 Hex | Y | | 4.1 | 2.5 Hex | Y |
| | 4.5 | 2.5 Hex | | 4.5 | 2.5 Hex | Y | | 4.3 | 2.5 Hex | Y |
| | 5.0 | 2.5 Hex | | 5.0 | 2.5 Hex | Y | | 4.8 | 2.5 Hex | Y |
| | 6.0 | 2.5 Hex | | 6.0 | 2.5 Hex | Y | | | | |
| INNO Int. | 3.7 | 3.1 Octa | Any Ridge Implant (sub type) | 3.5 | 2.3 Hex | N | i-clean | 3.7 | 3.1 Octa | Y |
| | 4.0 | 3.1 Octa | | 4.0 | 2.3 Hex | N | | 4.1 | 3.1 Octa | Y |
| | 4.5 | 3.1 Octa | | 4.5 | 2.3 Hex | N | | 4.3 | 3.1 Octa | Y |
| | 5.0 | 3.1 Octa | | 5.0 | 2.3 Hex | N | | 4.8 | 3.1 Octa | Y |
| | 6.0 | 3.1 Octa | | 5.5~8.0 | 2.3 Hex | N | | 6.0 | 3.1 Octa | Y |
| INNO Ext. | 3.7 | 2.7 Hex | | | | | e-clean | 3.5 | 2.4 Hex | N |
| | 4.0 | 2.7 Hex | | | | 4.1 | | 2.7 Hex | Y | |
| | 4.5 | 2.7 Hex | | | | 5.1 | | 3.4 Hex | Y | |
| | 5.0 | 3.4 Hex | | | | | | | | |
| | 6.0 | 3.4 Hex | | | | | | | | |